


Daily Doser

Automatic Pill Dispenser


The Daily Doser allows a doctor to take control of their patients medications. It is meant to help the people who may get confused with keeping track of taking multiple pills at different times of the day.



Operated by a battery, the doctor can program each dosage into the machine and have the specific pill box unlock at the time selected.

1. The doctor or pharmacist must insert their ID card into the side of the machine to gain access to the patient's settings and open all boxes.



They can easily program each pill release and alarms using the buttons underneath the display monitor.




- Left Moves mouse to the left.
- Right Moves mouse to the right.
- Menu Pulls up the menu for the settings options (Date/Month/Time).
- Enter Enter/End button.

2. Patients cannot unlock the boxes after the doctor has set the pills.

The patient must listen for the alarm to sound when it is time to take their pills.





- Alarm The alarm icon will turn red when the sound goes off.
- Enter Press enter to turn off the alarm.



The boxes set for that time will unlock, making it easy for the patient to know which pills to take at that time.

3. Close the door to the open pill box. Make sure the light turns red so that it is locked.



4. Take the pill box back to the pharmacist for refills.

